

Diagnostikleuchten | Diagnostic lights

FÜR DIE PRAXIS | FOR THE CLINIC

LED



Für Augenuntersuchungen
geeignet nach DIN EN 62471
suitable for eye examinations
in accordance with DIN EN
62471 standards



Cliplight Diagnostikleuchte LED | Cliplight diagnostic light LED

- Für Augenuntersuchungen geeignet nach DIN EN 62471
- Metallgehäuse, lackiert
- Batterien im Lieferumfang enthalten
- in praktischer Blisterpackung
- Betriebsspannung: 3V
- Batterielebensdauer ca. 50 Stunden
- Nickelfrei
- Beleuchtungsstärke ca. 24.000 Lux
- Lebensdauer LED ca. 100.000 Stunden
- Länge ca. 13 cm
- VE = 6 St.

- suitable for eye examinations in accordance with DIN EN 62471 standards
- made of metal, painted
- includes replaceable batteries
- practically packaged in blister pack
- operating voltage: 3V
- battery life: approx. 50 hours
- nickel-free
- light intensity approx. 24.000 Lux
- LED life approx. 100.000 hours
- length: approx. 13 cm
- PU = 6 items

rot red		REF 12.05401.014	lila violet		REF 12.05401.094
schwarz black		REF 12.05401.024	blau blue		REF 12.05401.234
gelb yellow		REF 12.05401.054	grau grey		REF 12.05401.264

VE = 6 St.

PU = 6 items

farblich sortiert | assorted colours  REF 12.05411.004

Diagnostikleuchte, weiß Diagnostic light, white

- aus Metall
- mit Druckknopf
- inkl. auswechselbaren Batterien 1,5 V Micro (Typ AAA)

- made of metal
- with push button
- includes replaceable batteries 1.5 V (type AAA)





VE = 5 St. | PU = 5 items

weiß | white

REF 12.05100.071

KaWe DIALIGHT XL Diagnostikleuchte KaWe DIALIGHT XL diagnostic light

- Kunststoffkopf
- mit Xenon-Lampe 2,5 V (REF 12.75145.003)
- Batterien 1,5 V, Mignon (Typ AA) 
- in Tasche

- plastic head
- with xenon bulb 2.5 V (REF 12.75145.003)
- batteries 1.5 V, mignon (type AA) 
- with pouch



Länge ca. 13 cm, Ø ca. 2 cm | length approx. 13 cm, Ø approx. 2 cm

night  REF 12.05302.021

ZUBEHÖR | ACCESSORIES

Spatelhalter | Tongue depressor holder

Spatelhalter zum einfachen Befestigen von Mundspateln an den Cliplight Diagnostikleuchten und Diagnostik-Leuchte, weiß.

Tongue depressor holder for easy attachment of wooden blades to the Cliplight diagnostic lights and diagnostic light, white.



VE = 5 St. | PU = 5 items

REF 12.05290.001

Tasche, allein | Pouch, separate

für KaWe DIALIGHT XL Diagnostikleuchten
REF 12.05302.021, 12.05302.231

for KaWe DIALIGHT XL diagnostic lights
REF 12.05302.021, 12.05302.231

schwarz | black

REF 01.78401.021

