



DIATROM®160

Dim.: 25,4 x 28,8 x 10,4 h cm

EN Electro-surgical unit

Monopolar and bipolar radio-surgical unit suitable for micro-coagulative operations, small and medium surgery. Different functions can be selected from the control panel: CUT, BLEND, COAG FORCED, COAG SOFT, BIPOLAR. Its digital display of the delivered power guarantees reliability and simple use.

FR Bistouri électrique

Bistouri électrique monopolaire/bipolaire conçu pour les interventions de micro coagulation, de petite et moyenne chirurgie. Son panneau de commande permet de sélectionner la fonction désirée: CUT, BLEND, COAG FORCED, COAG SOFT, BIPOLAR. L'affichage numérique de la puissance délivrée rend les actes sûrs et reproductibles.

IT Elettrobisturi

Elettrobisturi monopolare / bipolare progettato per micro coagulazione, piccola e media chirurgia. Il pannello di controllo consente di selezionare semplicemente la funzione desiderata: CUT, BLEND, COAG FORCED, COAG SOFT, BIPOLAR. Il display digitale della potenza erogata rende gli atti sicuri e riproducibili.

Technical specifications / Spécifications techniques / Specifiche tecniche

CODE	10100.30 NEW
Frequency / Fréquence / Frequenza	600 kHz
Max power / Puissance max / Potenza max.	CUT 160 W - 250 Ohm
	BLEND 120 W - 200 Ohm
	FORCED COAG 100 W - 150 Ohm
	SOFT COAG 80 W - 100 Ohm
	BIPOLAR 60 W - 100 Ohm
Patient circuit / Électrode neutre / Elettrodo neutro	F
Power supply / Alimentation / Alimentazione	115-230Vca / 50-60Hz
Power max absorption / Puissance max. absorbée / Potenza assorbita	350 VA
Weight / Poids / Peso	5 kg

Standard accessories / accessoires / accessori

		Qty
00100.03	Cord, 2 m / Câble 2 m / Cavo 2 m	1
F4243	Reusable sterile handpiece / Porte-électrodes monopolaire stérile / Portaelettrodo monopolare sterile	1
00404.08	Cable for cylindrical electrode / Câble de connexion électrode / Cavo di collegamento dell'elettrodo	1
5365	Small neutral plate / Plaque de référence / Piastra di riferimento	1
00304.00	Pedal switch / Pédale interrupteur / interruttore a pedale impermeabile	1
00500.00	Kit of 10 electrodes / Kit de 10 électrodes assorties / Kit di 10 elettrodi assortiti	1
MA280	User's manual / Notice d'emploi / Manuale d'uso	1