

# precisa<sup>®</sup> N

Precision that lasts a lifetime



**↓ BIHS**

Handheld sphygmomanometer with sturdy aluminium or ABS enclosure. Available with the choice of single or twin tube bladder cuffs. Also optionally available in a latex-free version (when ordering enter LF for latex-free).

- Aluminium/Metal version with aluminium enclosure and tube plug-in connector
- Precision air release valve with fine adjustment and wear-free.
- Microfilter protects air release valve and measuring system.
- Specially tempered copper-beryllium diaphragm, with precision movement and almost non-aging.
- Diaphragm pressure-loading capacity up to 600 mm Hg.
- Metal valve base.
- Latex bulb for pumping up the cuff.
- Ergonomic spoon-shaped grip of stainless steel.
- Aluminium casing, stainless steel spoon.
- Linear, tilted aluminum scale, from 0 to 300 mm Hg (easily legible).
- Maximum error tolerance +/- 3 mm Hg.
- Comprehensive selection of different sizes and versions of cuffs.

<b>No. 1360-101</b>	precisa <sup>®</sup> N, plastic/metal, 1-tube, bandage cuff (adult)	<b>No. 1361-103</b>	precisa <sup>®</sup> N, double tube, bandage cuff (adult)
<b>No. 1360-104</b>	precisa <sup>®</sup> N, plastic/metal, 1-tube, hook cuff (adult)	<b>No. 1361-106</b>	precisa <sup>®</sup> N, double tube, hook cuff (adult)
<b>No. 1360-107</b>	precisa <sup>®</sup> N, plastic/metal, 1-tube, velcro cuff (adult)	<b>No. 1361-109</b>	precisa <sup>®</sup> N, double tube, velcro cuff (adult)
<b>No. 1362-101</b>	precisa <sup>®</sup> N, aluminium, 1-tube, bandage cuff (adult)	<b>No. 1363-103</b>	precisa <sup>®</sup> N, aluminium, double tube, bandage cuff (adult)
<b>No. 1362-104</b>	precisa <sup>®</sup> N, aluminium, 1-tube, hook cuff (adult)	<b>No. 1363-106</b>	precisa <sup>®</sup> N, aluminium, double tube, hook cuff (adult)
<b>No. 1362</b>	precisa <sup>®</sup> N, aluminium, 1-tube, velcro cuff (adult)	<b>No. 1363</b>	precisa <sup>®</sup> N, aluminium, double tube, velcro cuff (adult)

\* Further items are available, please contact your Riester distributor for further information.

Please visit **riester.de** for more information