

# Lightweight rollator Gigo 2G

Item No. 735100200



- Back rest now removable and foldable
- Particularly light and attractive aluminum rollator
- One-hand folding mechanism ensures safe and comfortable folding
- Ergonomic handles can be adjusted to 6 height positions
- Max. load of basket and tray 5 kg each
- Supplied as standard with basket, tray, backrest, cane holder and tilt aid

## Technical Data

Order number	735100200	Width between handles	40,5 cm
Colour	silver	Folded width	61,5 cm
Overall width	61,5 cm	Folded height	100 cm
Overall length	61 cm	Folded length	28 cm
Handle height	87 - 96 cm	Total weight	7.8 kg
Seat width	34 cm	Weight without accessories	7 kg
Seat depth	33.5 cm	Max. Load	120 kg
Seat height	- 54 cm	PU Pal.	10 Pal.

# Lightweight rollator Gigo 2G

Item No. 735100200

**drive**

**DeVilbiss**  
HEALTHCARE

## Accessories



**Basket**

Item No.710000001



**Tray**

Item No.710000002



**Rollator bag**

Item No.710000075



**Rollator bag**

Item No.710000076



**One-hand brake**

Item No.710000081



**Cane holder Gigo / Gigo 2G**

Item No.710000090