

Easy READER+®

New Generation Quantitative Test Reader

The VEDALAB Easy Reader+® POC reader allows quantitative measurements from fingerprick whole blood samples using lateral flow tests.

- Flexible user interface SW :
Easy mode and advanced mode
- Stand alone HW :
 - Integrated printer
 - Touch screen
 - Barcode reader (optional)
 - CCD camera
- Optimal communication :
 - Interfaces : USB A, USB B, Ethernet, RS232
 - Supported protocols : HL7, LIS2-A2
- Data storage : last 1,000 results with images
- Adaptable system : Multilingual, web upgrade SW, usable with AA batteries (portable)
- Cost effective : affordable instrument usable with a very wide range of tests... and many more to come !
- Many languages available (English, French, German, Spanish, etc.)
- Possibility of qualitative tests measurements

Available quantitative tests

- | | | | |
|--|------------------|---|---------------|
| ■ AFP | ■ Cystatin C | ■ HS-Insulin | ■ Transferrin |
| ■ CA-125 | ■ C-Peptide | ■ LH | ■ Troponin I |
| ■ CA 15-3 | ■ D-Dimer | ■ Microalbumin | ■ TSH |
| ■ CA 19-9 | ■ FABP | ■ Myoglobin | ■ HS-TSH |
| ■ Calprotectin | ■ Ferritin | ■ NT-proBNP New! | ■ ULTRA-TSH |
| ■ CEA | ■ FOB | ■ Procalcitonin | |
| ■ CK-MB | ■ FSH | ■ Prolactin | |
| ■ Cortisol | ■ HbA1c | ■ HS-Prolactin | |
| ■ COVID-19 spike IgM/IgG New! | ■ Hb-Total | ■ PSA (quantitative or semi-quantitative) | |
| ■ CRP (laboratory or doctor's office) | ■ hCG (S/P or U) | ■ PSA Free | |
| ■ HS-CRP (laboratory or doctor's office) | ■ beta-hCG | ■ Total T3 | |
| | ■ IgE | ■ Total T4 | |
| | ■ Insulin | | |

Coming soon

- Progesterone
- Testosterone
- Vitamin D

Available qualitative tests

- | | | | |
|----------------|------------------------------------|-----------------------------------|---------------------|
| ■ Adenovirus | ■ C.d. Combo Toxin A+B New! | ■ Combo Influenza A+B New! | ■ Strep. A |
| ■ BNP | ■ hCG (10 min) | ■ Legionella | ■ Strep. Pneumoniae |
| ■ C.d. GDH | ■ hCG (5 min) | ■ Rotavirus | ■ Tetanus |
| ■ C.d. Toxin A | ■ Influenza A | ■ Combo Rota-Adeno New! | ■ Tuberculosis |
| ■ C.d. Toxin B | ■ Influenza B | ■ RSV | |

■ Most of items are CE certified



VEDA-LAB

Rue de l'Expansion, ZAT du Londeau - Cerisé - B.P. 181
61006 Alençon Cedex France
Tél. : (33)2 33 27 56 25 - Fax : (33)2 33 27 70 60
E-mail : veda.lab@wanadoo.fr
www.vedalab.com

