

MIO-SONIC

Mio-Sonic is the compact and powerful ultrasound device, very suitable for muscle contractures and joint stiffness.

MEDICAL DEVICE CLASS IIa

MIO-SONIC is a medical device CE0197



KIT

- ✓ Mio-Sonic
- ✓ Medical power supply
- ✓ Ultrasound gel
- ✓ User manual
- ✓ Carriage and storage bag
- ✓ Box

TECHNICAL FEATURES

- ✓ Professional device for ultrasoundtherapy
- ✓ 1 MHz ultrasound, maximum output power 6,4W, maximum intensity 1,6W/cm²
- ✓ Collimated ultrasound beam for an optimal penetration beam
- ✓ 3 ultrasound beam intensity levels: 0,08 W/cm² (low) – 0,8 W/cm² (medium) – 1,6 W/cm² (high)
- ✓ Medical power supply
- ✓ Direct contact device (not suitable for water immersion)

DIMENSIONS 200 x 47 x 69 mm

WEIGHT 190 g

MIO-SONIC

PROGRAMS

TREATMENTS

Osteoporosis	Tendinitis	Bursitis
Arthritis	Pubalgy	Brachialgy
Neuralgias	Tallonitis	Brachialgy
Cervicalgias	Contusions	Chondropathy
Sciatalgy	Sprains	Cruralgy
Lumbalgy	Hematoma	Cramps
Carpal Tunnel	Oedema	Venous stagnation
Epycondilitis	Lymphedema	Cramps
Epitrochleitis	Adhesions	
Periarthritis	Algodystrophy	

HIGH PERFORMANCE DEVICE

Innovative and high performance device both for home use and professional use. It combines compactness and high performances in a small, light and easy to use applicator. The lightness and handiness allows a continuous use without any difficulty.

Mio Sonic can be used for the following applications:

- ✓ Pain therapy
- ✓ Fitness treatments

MIO-SONIC

MANUFACTURER

Shenzhen Dongdixin Technology Co., Ltd.
Floor 1-2, No.3 Building, Fanshen Xusheng
Industrial Estate Xilixiaobaimang 518108
Nanshan District, Shenzhen P.R. China

EC REPRESENTATIVE

Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany

DISTRIBUTED BY

I.A.C.E.R Srl Via Enzo Ferrari 2 – 30037 Scorzè (VE) Italy