



AUDITEST



The manual audiometer **AudiTest** is achieving an **accurate screening** of hearing faculties. Convenient and performant, its size and ergonomics allow the Healthcare professional to use it with just one hand.



AUDITEST



Lightweight, portable, **manageable**. It can be used with just one hand.



Perfect **ergonomics** thanks to its quiet keys and its LCD display



- **Continuous or pulsed** sound
- **Frequencies** from 125 to 8000 Hz, output level range from -10 to 100 HLdB



Perfectly adapted to the **screening in schools**



- **Device designed and manufactured in France**
- **Warranty : 3 Years ***



End user manual : English, German, Italian, Spanish, Portuguese, Dutch, Hungarian, Polish, Romanian, French

WARRANTY 3 YEARS (EXCEPT FOR THE HEADPHONES) – FRENCH DESIGN AND MANUFACTURING

TECHNICAL INFORMATION

- **Type** : 4 – Screening and Control (according to EN 60645-1) Medical Marking EC 0459 – Class IIa
- **Output** : AC – Air Conduction
- **Frequencies range (11 values)** : 125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000 Hz (+/- 1%)
- **Sound level** : from -10 HLdB to 100 HLdB

- **Increment of sound level** : 5 dB
- **Power supply** : 2 battery sticks 1,5V -Type AA / LR06 – Supplied
- **Stimulus** : continuous or pulsed sound
- **Headphones** : light, comfortable
- **Dimensions of the audiometric box** : 79 x 122 x 33 mm
- **Packaging** : soft kit
- **Total weight of the kit** : 600 g (155 g for the audiometric box only).

Electronica Medical offers a complete range of devices for screening the hearing losses, which are dedicated to the educational and occupational medicine.

Electronica Technologies reserves itself the right to modify the specifications indicated without prior notice. The photos are not contractual.

Headquarters and Industrial Center : Chemin du Pré Long – 03200 Abrest – France. Tél. +33 (0)4 70 96 52 90 - Fax +33 (0)4 70 96 01 23 - contact@electronica-medical.com

* Except for the headphones