



# Height rod, foldable

5003

- ▶ Ultrasound measuring technology, temperature-compensated
- ▶ Tilt sensor for correction if measuring instrument is not horizontal, resulting in an error message
- ▶ Low space requirement due to fold function
- ▶ Units in cm or inch



## Features

The length measuring rod from Soehnle Professional, which is particularly suitable for paediatricians and hospitals, not only impresses with its simple, stylish design, but also with an integrated automatic electronic spirit level. This enables free measurement in space and guarantees the highest measuring precision and speed at all times. Another advantage of the length measuring rod is its unique folding function, which means it takes up hardly any space and is easy to stow away. Moreover, thanks to an acoustic acknowledgement and automatic fixing of the values after the measurement has been completed, the process does not need to be laboriously monitored.

- ▶ Ultrasound measuring technology, temperature-compensated
- ▶ Tilt sensor for correction if measuring instrument is not horizontal, resulting in an error message
- ▶ Hold function for easy reading of the height value after measurement
- ▶ Low space requirement due to fold function
- ▶ Plastic housing
- ▶ LCD display, digit height 9 mm
- ▶ Acoustic acknowledgement at end of measurement and reading displayed
- ▶ Units in cm or inch
- ▶ Battery operation with 2 x 1,5 V AA batteries (included in scope of delivery)

## Products

Order number Variation	Approvable	Division mm	Measuring range cm	Dimensions (mm)
<b>5003.01.001</b> Height rod, foldable	No	10	50 - 240	345×92×28mm