



Baby scale "Multina Comfort"

8352

- ▶ Simple operation with just one button
- ▶ With removable baby tray
- ▶ Can be used for children up to 20 kg
- ▶ Analogue/digital baby length measuring rod optional (see accessories)



Features

The Soehnle Professional Baby Scale 8352 convinces by its simple but high-class design in combination with highest quality - made in Germany. The intuitive one-button operation, the tare function and the easy-to-clean construction make the baby scale a popular companion for midwives, mothers and birth centres.

The highlight of this scale is the removable baby car seat. This means that even small children up to 20 kg can be weighed on the scale in a standing or sitting position.

- ▶ Simple operation with just one button
- ▶ With removable baby tray
- ▶ Easy-to-read LC display, 13 mm digit height
- ▶ Tare function: the weight can be tared by placing a weight on it
- ▶ Hold function for convenient reading of the weight value after weighing
- ▶ Can be used for children up to 20 kg
- ▶ Unit changeover switch kg/lb
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 9V block (included in scope of delivery)
- ▶ Battery status display
- ▶ Analogue/digital baby length measuring rod optional (see accessories)

Products

Order number Variation	Approvable	Division g	Maximum load kg
8352.01.001 Baby scale "Multina Comfort"	No	10	20

Accessories

Order number

Variation

Description

5005.11.001



Baby height rod

5040.03.001



Carry bag black, with detachable shoulder strap and durable zipper

5020.01.001

Baby tray (replacement)